

Article Freka MCT Oil Ampules

Art. No. FK 3087364

Intended use: Lubricant for removing guide wires

Class: I

Rule: 1

Number of pages 02

Date of revision 12.07.2011

Replaces revision from


Issued by A. Manck
Function Manager Regulatory Affairs

Signed


.....

Revision list

Date	No. of revision	Modification
12.07.2011	0	First release



Article Freka MCT Oil Ampules

Art. No. FK 3087364

Parts List

Position No.	Description	Material
1	MCT Oil 3ml	Caprylic- / capric triglyceride
2	Ampule	PE [Polyethylene
Packaging		
	PE-Beutel / PE-bag	Polyethylene
	Karton / carton	Cardboard

Drawing

