


Article Number 4940001

PharmaNeedle

Description Capped side port needle with luerlock connection (Sterilized)



Length: 76 mm (incl. cap)
 46 mm (without cap)
Needle Ø: 1.5 mm / 17G
Needle shaft: 20 mm
Needle material: Stainless steel 304
Packaging 10 pcs in Soft blister
 25 blisters in Inner Carton
 18 Cartons in Outer Box

Product label


PharmaNeedle. REF XXXXXX
 LOT XXYXX-X
 Ø 1,5/ 17G | 20 mm
 Not for human use 10 pcs.
 XXXX-XX
 XXXX-XX
 CE0123
 Süddeutsche Feinmechanik GmbH
 Brückenstraße 5
 D-63607 Wächtersbach
 Medical Dispensing Systems
 STERILE EO
 2

Inner Carton Label


Inner carton dimensions: 27,5 x 14 x 12,5 cm



PharmaNeedle. REF XXXXXX
 LOT XXYXX-X
 Ø 1,5/ 17G | 20 mm
 Not for human use 250 pcs.
 XXXX-XX
 XXXX-XX
 CE0123
 Süddeutsche Feinmechanik GmbH
 Brückenstraße 5
 D-63607 Wächtersbach
 Medical Dispensing Systems
 STERILE EO
 2

Outer Box Label

Outer box dimensions: 59 x 45 x 39,5 cm



PharmaNeedle. REF XXXXXX
 LOT XXYXX-X
 Ø 1,5/ 17G | 20 mm
 Not for human use 4500 pcs.
 XXXX-XX
 XXXX-XX
 CE0123
 Süddeutsche Feinmechanik GmbH
 Brückenstraße 5
 D-63607 Wächtersbach
 Medical Dispensing Systems
 STERILE EO
 2

Contraindications

PharmaNeedle is to be used in combination with PharmaHelp only

The PharmaNeedle is not for human use

Sterilization Method

Ethylen Oxide

Shelf life

5 years from the month of production

Certificate

CE