

Article Number	4945101	PharmaHelpBag
Description	Sterilized waste bag to collect surplus carrier solution in PharmaHelp	
	Nominal capacity:	300 ml
	Material:	Multi-layer polyolefin film
	Packaging:	Single packed in sealed soft blister
		75 soft blisters per box
	Box dimensions 29 x 27 x 24 cm	
Product label		
Box Label		
Contraindications	Bag is not intended to carry medication to be delivered to a patient The PharmaHelpBag is not CE-Certified	
Sterilization Method	Y—Irradiation using VDmax method. A Typical dose is 25Kgy	
Shelf life	5 years from the month of production	
Certificate	Not CE-Certified	